

# Laboratory Approach to Haematological Disorders

Dr. Albert Sin  
Clinical Assistant Professor  
Department of Pathology  
Faculty of Medicine  
The University of Hong Kong





# Learning Objectives


- Understand the scope of haematological disorders
- Describe laboratory tests useful for investigation
- Rationalise the use of various tests in patient management



# Elements of the Haematological System



## ❑ Organs

- Bone marrow
- Thymus
- Spleen
- Lymph nodes




## ❑ Cells

- White blood cells
- Red blood cells
- Platelets



## ❑ Plasma

- Coagulation factors
- Anticoagulants






# Elements of the Haematological System

## ❑ Organs

- Bone marrow
- Thymus
- Spleen
- Lymph nodes




## ❑ Cells

- White blood cells
- Red blood cells
- Platelets



## ❑ Plasma

- Coagulation factors
- Anticoagulants





# Nature of Haematological Disorders

- Hereditary vs acquired
- Benign/reactive vs malignant
- Primary vs secondary (to systemic diseases)



# Laboratory Investigations for Haematological Diseases

- Aim of clinical laboratory investigation
  - Establish diagnosis, predict prognosis, guide treatment, help monitoring



Collect Date :	05/04/03
Collect Time :	--:--
Receive Date :	05/04/03
Receive Time :	18:00
Request No. :	HN021608
Urgency :	URGENT
Complete Blood Count	
WBC	8.7
RBC	2.96 L
HGB	8.3 L
HCT	0.248 L
MCV	83.8
MCH	27.9
MCHC	33.3
RDW	13.5
PLT	218
MPV	9.8
Film Review	N



# Basic Principles of Clinical Laboratory Investigation



- Simple fast screening to complex confirmatory tests
- Cheap to expensive tests
- Non-invasive to invasive tests



# Blood Cell & Bone Marrow Investigations



# Complete Blood Counting - The Simplest Haematology Test

- Simple, fast, cheap, non-invasive and excellent screening
- Indicate cell count abnormalities and average red cell characteristics
- Non-specific as to the underlying diagnosis



Collect Date :	05/04/03
Collect Time :	--:--
Receive Date :	05/04/03
Receive Time :	18:00
Request No. :	HN021608
Urgency :	URGENT
-----	
Complete Blood Count	
WBC	8.7
RBC	2.96 L
HGB	8.3 L
HCT	0.248 L
MCV	83.8
MCH	27.9
MCHC	33.3
RDW	13.5
PLT	218
MPV	9.8
Film Review	N



# Automated Blood Counting





# What is the Cause of Anaemia?

- Get the most information from simple screening tests - Complete Blood Counting

- Red cell indices
- White cells and platelets
- Reticulocyte count
- Peripheral blood smear

Collect Date :	05/04/03
Collect Time :	--:--
Receive Date :	05/04/03
Receive Time :	18:00
Request No. :	HN021608
Urgency :	URGENT

## Complete Blood Count

WBC	8.7
RBC	2.96 L
HGB	8.3 L
HCT	0.248 L
MCV	83.8
MCH	27.9
MCHC	33.3
RDW	13.5
PLT	218
MPV	9.8
Film Review	N



# CBP Report

Collect Date : 05/04/03  
Collect Time : --:--  
Receive Date : 05/04/03  
Receive Time : 18:00  
Request No. : HN021608  
Urgency : URGENT

## Complete Blood Count

WBC	8.7
RBC	2.96
HGB	8.3
HCT	0.248
MCV	83.8
MCH	27.9
MCHC	33.3
RDW	13.5
PLT	218
MPV	9.8
Film Review	N





# Clinical Classification of Anaemia

## ■ By red cell size (Mean Cell Volume)

### Microcytic

Iron deficiency

Thalassaemia

### Normocytic

Haemolysis

Anaemia of chronic  
disease

Renal failure

### Macrocytic

Megaloblastic anaemia

Aplastic anaemia

Myelodysplasia


Liver disease



# Peripheral Blood Smear Examination

## - A Useful Supplement to CBC

- Confirm blood counts, assess morphology
- Provisional diagnosis
- Guide further special investigation




Collect Date :	05/04/03
Collect Time :	--:--
Receive Date :	05/04/03
Receive Time :	18:00
Request No. :	HN021608
Urgency :	URGENT
-----	
Complete Blood Count	
WBC	8.7
RBC	2.96 L
HGB	8.3 L
HCT	0.248 L
MCV	83.8
MCH	27.9
MCHC	33.3
RDW	13.5
PLT	218
MPV	9.8
Film Review	N



# Blood Smear in Iron Deficiency Anaemia

Small and pale red cell

Pencil cell





# Confirmation of Iron Deficiency

## Biochemical findings


- low serum iron (usable iron)
- high serum transferrin (compensatory response aiming to increase transport)
- low transferrin saturation
- low serum ferritin (body iron storage)



# Blood Smear in Thalassaemia Trait

Small red cell


Target cell





# $\alpha$ -Thalassaemia Diagnosis


## Detection of HbH ( $\beta_4$ ) Inclusion Bodies






# $\beta$ -Thalassaemia Diagnosis

## High Performance Liquid Chromatography to Measure Hb A<sub>2</sub>






# Some peripheral blood smear examination results suggest a bone marrow disease



circulating blasts





unexplained cytopenia



# Bone Marrow Examination

- Invasive but safe, requires expertise in assessment








# Bone Marrow Examination

- Determine whether the marrow is the primary site of abnormality
  - Inadequate vs ineffective haemopoiesis in cytopenia

- Presence of abnormal cells - blasts, infiltrates, dysplasia



# Bone Marrow Examination






congenital dyserythropoietic anaemia

aplastic anaemia

miliary tuberculosis



# Bone Marrow Examination



acute promyelocytic leukaemia


myelodysplastic syndrome

primary myelofibrosis



# Bone Marrow Examination

- Provide optimal material for further special investigations in malignant cases



# Special Investigations for Malignant Haematological Diseases




- To classify malignant disease for choice of specific treatment
  - Different myeloid types, B, T, natural killer cell
- To predict prognosis for counselling and risk-adapted therapy
- To monitor residual disease after treatment
  - Minimal residual disease (MRD) not detectable by morphology

# Special Investigations for Malignant Haematological Diseases



## ■ Cytochemistry

- Classification of acute leukaemia
- Demonstrate enzymatic function of blasts




# Special Investigations for Malignant Haematological Diseases



## ■ Immunophenotyping

- Classification of leukaemias
- Using antibodies to detect lineage-associated antigens in abnormal cells





Haemopoietic precursors: CD34, Tdt  
B-lineage: CD19, CD20, CD22, CD79a  
T-lineage: CD3, CD4, CD5, CD8, CD7, TCR  
Myeloid: CD13, CD33, MPO  
Megakaryocytic: CD41, CD42b, CD61  
Erythroid: glycophorin A

# Special Investigations for Malignant Haematological Diseases



- Immunophenotyping
  - Flow cytometer - automated machine for fast and detailed immunophenotyping





# Special Investigations for Malignant Haematological Diseases



## ■ Immunophenotyping



- Laser light excites fluorochrome on antibodies bound to abnormal cells



# Special Investigations for Malignant Haematological Diseases



- Immunophenotyping
  - MRD monitoring if leukaemic cells have a specific immunophenotype





# Genetic Testing for Haematological Diseases





- As confirmation for inherited haematological diseases with known genetic defects (e.g. thalassaemia)
- Malignant haematological diseases are also genetic diseases
  - Acquired mutations

# Genetic Testing for Malignant Haematological Diseases



## ■ Cytogenetic studies

- Diagnosis and classification of leukaemias
- Culture malignant cells and examine their chromosomes





Philadelphia chromosome

# Genetic Testing for Malignant Haematological Diseases



## ■ Cytogenetic studies

- Prediction of prognosis
- Culture malignant cells and examine their chromosomes






# Genetic Testing for Malignant Haematological Diseases



## ■ Molecular genetic tests


- Classification, prediction of prognosis and MRD monitoring
- Use sensitive molecular techniques to determine lineage of abnormal cells, gene mutation and gene fusion





# Genetic Testing for Malignant Haematological Diseases

- Molecular genetic tests for MRD monitoring




## Clinical Suspicion

## Blood Counts & Smear Morphology

Collect Date : 05/04/03  
Collect Time : 11:45 AM  
Receive Date : 05/04/03  
Receive Time : 18:00  
Request No. : HEMO 1059  
Urgency : URGENT


Complete Blood Count

Item	Value
WBC	8.7
RBC	2.96 L
HGB	9.3 L
HCT	0.288 L
MCV	83.8
MCH	27.9
MCHC	31.3
RDW	13.5
PLT	218
NWBC	0.8
#/ml Review	N



## Blood Tests

## Bone Marrow Morphology



## Non-malignant Dx

## Cytochemistry / Immunophenotyping




## Cytogenetics / Molecular Tests



## Malignant Dx

## MRD Detection



## Immunophenotyping / Molecular Tests

treatment



# Bleeding & Thrombosis Investigations



# You Have a Patient with Bleeding






# Components of Haemostasis

- Blood vessel
- Platelets
- Coagulation factors
- von Willebrand factor (VWF)
- Natural anticoagulants (protein C, protein S, antithrombin)

Blood clot formation






# Coagulation Factors, VWF & PC, PS AT Are in Patient's Plasma

citrate bottle



anticoagulated  
blood






# Coagulometers





# Prothrombin Time (PT)



Abnormal PT

factor VII

factor X

factor V

prothrombin

fibrinogen



# Prothrombin Time

anticoagulated patient plasma

+

thromboplastin (tissue factor & phospholipid)

+


calcium



= PT (sec)



# Activated Partial Thromboplastin Time (APTT)



## Abnormal APTT

factor XII

factor XI

factor IX

factor VIII

factor X

factor V

prothrombin

fibrinogen



# Activated Partial Thromboplastin Time

anticoagulated patient plasma

+

contact activator

+


phospholipid & calcium



= APTT (sec)



# Thrombin Time (TT)



Abnormal TT

fibrinogen



# Thrombin Time

anticoagulated patient plasma

+

thrombin

+



= TT (sec)



# Other Tests in Blood & Plasma for Bleeding and Thrombotic Diseases

- Depending on the clinical presentation and screening test results
  - Specific coagulation factor assay
  - Coagulation inhibitor testing
  - VWF testing
  - PC, PS and AT assay
  - Platelet function test



# Conclusions

- Laboratory investigations provide disease diagnosis, prognosis and monitoring
- Rational use of special haematology tests allow cost-effective use of resources for patient management



# Photo Credits and References

- Interactive Haematology Image Bank. B. Bain. Blackwell Science.
- Wintrobe's Atlas of Clinical Hematology. D. Tkachuk & J. Hirschmann. Lippincott Williams & Wilkins.
- Color Atlas of Clinical Hematology. 4<sup>th</sup> Ed. V. Hoffbrand, J. Pettit & P. Vyas. Elsevier.
- Diagnostic Techniques in Hematological Malignancies. W. Erber. Cambridge University Press.
- Essential Haematology by AV Hoffbrand, PAH Moss, and JE Pettit 5th Edition 2006, Blackwell Publishing